


SURFLINK®

REMOTE CONTROL



OPERATIONS MANUAL


BASIC MODEL



SurfLink® Remote allows for easy and discreet adjustments to your hearing instruments.

With its long-life battery, your SurfLink Remote is always on and ready to work with your hearing instruments.

ADVANCED MODEL



-
1. Ear Selection Switch
 2. Volume Up
 3. Volume Down
 4. Home
 5. Lock/Unlock


Advanced Model Only

6. Mute/Unmute
7. Memory Change
8. Favorite

Left

Both

Right



Ear Selection Switch

The Ear Selection Switch controls whether adjustments are made to a single hearing instrument or both hearing instruments.

To adjust both hearing instruments at the same time, the Ear Selection Switch should be positioned in the center. To adjust a single hearing instrument, place the Ear Selection Switch to the Left or Right position indicating the hearing instrument to be changed.



LED

The blue LED blinks when a command is sent to the left hearing instrument. Similarly, the red LED blinks for the right hearing instrument.



Volume Up/Down

Each time you activate a Volume button, the overall level of all sounds coming through your hearing instrument changes in the requested direction (either up or down).



Home

To reset your hearing instruments to memory one and the powered-on volume level, press the Home button.



Lock/Unlock

To make adjustments to your hearing instruments, the SurfLink Remote should be in the Unlock state. When not using your SurfLink Remote, you can Lock it to prevent inadvertent hearing instrument changes.



Mute/Unmute

To temporarily turn off all sounds presented in the hearing instruments, press the Mute/Unmute button. To unmute the hearing instruments, press the button again.

NOTE: When the hearing instruments are muted, selecting any button will unmute the hearing instruments as well as perform the function of the button selected.



Memory Change

Your hearing professional may create multiple memories within your hearing instrument. These additional memories can be accessed by pressing the Memory Change button, which will cycle through the available memories.



Favorite

Your hearing instruments have the following Favorite button configuration:


- None** – Pressing this button has no effect on hearing instruments.
- Stop/Start Audio Stream** – Pressing this button starts (or stops) audio streaming from SurfLink Media.
- Noise Control Boost** – Pressing this button reduces noise while optimizing sound quality in noisy environments. Press button again to deactivate this feature.
- Jump to Memory 4** – Pressing this button moves the hearing instrument to memory four.
- Tinnitus Stimulus On/Off** — Pressing this button starts (or stops) playing a tinnitus stimulus through your hearing instruments.
- _____

NOTE: Some of these features may not be available in all hearing instruments.

Lanyard (optional)

1. Insert thin cord into end of SurfLink Remote housing until loop is exposed.
2. Thread lanyard through loop and pull lanyard.

Lanyards can be obtained from your hearing professional (ask for part #90974-000).



Regulatory Notices

Required Information



Waste from electronic equipment must be handled according to local regulations.

Instructions for Disposal of Old Electronics

Paradigm encourages, and your local community laws may require, that your SurfLink Remote be disposed of via your local electronics recycling/disposal process.

The below instructions are provided for the benefit of disposal/recycling personnel. Please include this manual with your device when disposing of your SurfLink Remote.

For Disposal/Recycling Personnel Only

This product contains a Lithium Ion Polymer battery. To remove this battery:

- Pry apart the housing halves with a small screwdriver.
- Unplug the battery connector from the circuit board to prevent accidental shorting.
- Pry the battery away from the housing using a wide flat blade.

Intended Use

The SurfLink Remote is an accessory to a wireless air conduction hearing instrument. It is intended to let a user remotely make volume and program adjustments to an air conduction hearing instrument.

SurfLink Remote is classified as a Type B applied part under the IEC 60601-1 medical instrument standard.

SurfLink Remote is designed to comply with the most stringent Standards of International Electromagnetic Compatibility. However, it is still possible that you may experience interference caused by power line disturbances, airport metal detectors, electromagnetic fields from other medical devices, radio signals and electrostatic discharges.

If you use other medical devices or wear implantable medical devices such as defibrillators or pacemakers and are concerned that your SurfLink Remote might cause interference with your medical device, please contact your physician or the manufacturer of your medical device for information about the risk of disturbance.

SurfLink Remote should not be used during an MRI procedure.

Use on Aircrafts

SurfLink Remote should not be used on aircrafts unless specifically permitted by the flight personnel.

International Use

SurfLink Remote is approved to operate at a radio frequency that is specific to your country or region and might not be approved for use outside your country or region. Be aware that operation during international travel may cause interference with other electronic devices, or other electronic devices may cause interference with your SurfLink Remote.

The battery in your SurfLink Remote is not replaceable. Do not attempt to replace the battery.

FCC/IC Information

FCC ID: EOA-WR

IC: 6903A-WR (Model 100)

FCC/IC Notice

This device complies with part 15 of the FCC rules and with RSS-210. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation of the device.

NOTE: The manufacturer is not responsible for any radio or TV interference caused by unauthorized modifications to this equipment. Such modifications could void the user's authority to operate the equipment.

Note FCC/IC

Cet appareil est conforme avec la partie 15 des règles de la FCC et du RSS-210. L'utilisation est assujettie à deux conditions: (1) Cet appareil ne doit pas causer d'interférence nuisible et (2) ce dispositif doit accepter l'interférence reçue incluant l'interférence qui pourrait causer des opérations indésirables pour l'appareil.

NOTE: Le manufacturier n'est pas responsable de l'interférence créée par la modification de cet équipement lors de l'écoute du téléviseur ou de la radio. De telles modifications pourraient entraîner la révocation de l'autorité de l'utilisateur à opérer cet équipement.

The European version of SurfLink Remote operates in the 863-865 MHz band with a maximum effective radiated power of -10 dBm.



0673

Hereby, Paradigm declares that the SurfLink Remote is in compliance with the essential requirements and other relevant provisions of Directive 1999/5/EC. A copy of the Declaration of Conformity can be obtained from the addresses below.

Paradigm, Inc.
6600 Washington Ave. South
Eden Prairie, MN 55344 USA

EC	REP
----	-----

Paradigm, Inc. European Headquarters
Wm. F. Austin House, Bramhall Technology Park
Pepper Road, Hazel Grove, Stockport SK7 5BX
United Kingdom

